



PATIENT COOLING

PRE- & IN-HOSPITAL PRACTICAL INSIGHTS

CARDIAC ARREST IN VIENNA

Patient data*

- 7.030 patients without signs of circulation (206/100,000/year)
- 1.448 resuscitation attempts by the Vienna Ambulance Service
- 361 patients with sustained ROSC (25%)
- 164 (11.3%) discharged
- 126 (8.7%) with good outcome / CPC 1-2

* * 2009-2010

Out of Hospital Cardiac Arrest in Vienna: Incidence and Outcome.
Nürnberg et al. Resuscitation 2013





VIENNA – A SHORT STORY OF COOLING

HEART - LUNG RESUSCITATION

I FIRST AID: OXYGENATE THE BRAIN IMMEDIATELY

IF UNCONSCIOUS 1 or 2 operators

Airway - TILT HEAD BACK

IF NOT BREATHING

Breathe - INFLATE LUNGS 3-5 TIMES, MAINTAIN HEAD TILT

MOUTH-TO-MOUTH, MOUTH-TO-NOSE, mouth-to-adjunct, bag-mask

- FEEL PULSE
- IF PRESENT - CONTINUE LUNG INFLATIONS
- IF ABSENT -

Circulate - COMPRESS HEART ONCE A SECOND. ALTERNATE 2-3 LUNG INFLATIONS WITH 15 STERNAL COMPRESSIONS UNTIL SPONTANEOUS PULSE RETURNS.

for physicians only

II START SPONTANEOUS CIRCULATION

Drugs - EPINEPHRINE: 1.0 mg (10 CC OF 1:1000) I.V. OR 0.5 mg INTRACARDIAC. REPEAT LARGER DOSE IF NECESSARY.

SODIUM BICARBONATE: APPROXIMATELY 3.75 G/50 CC (1/2 DOSE IN CHILDREN) I.V. REPEAT EVERY 5 MINUTES IF NECESSARY

E. K. G. - FIBRILLATION: EXTERNAL ELECTRIC DEFIBRILLATION REPEAT SHOCK EVERY 1-3 MINUTES UNTIL FIBRILLATION REVERSED

IF ASYSTOLE OR WEAK BEATS: EPINEPHRINE OR CALCIUM I.V.

Fluids - I.V. PLASMA, DEXTRAN, SALINE

Do not interrupt cardiac compressions and ventilation. Tracheal intubation only when necessary. AFTER RETURN OF SPONTANEOUS CIRCULATION USE VASOPRESSORS AS NEEDED. e.g. NOREPINEPHRINE (Levophed) I.V. DRIIP

III SUPPORT RECOVERY (physician specialist)


Gauge EVALUATE AND TREAT CAUSE OF ARREST

Hypothermia START WITHIN 30 MINUTES IF NO SIGN OF CNS RECOVERY

Intensive Care SUPPORT VENTILATION: TRACHEOTOMY, PROLONGED CONTROLLED VENTILATION, GASTRIC TUBE AS NECESSARY

SUPPORT CIRCULATION
CONTROL CONVULSIONS
MONITOR

Figure 1. Heart-lung resuscitation (cardiopulmonary-cerebral resuscitation). First composition in 1961, Pittsburgh, PA. Reproduced with permission from Safar P. Community-wide CPR. J Iowa Medical Society 1964 (Nov); pp 629-635.



VIENNA RESUSCITATION APPROACH

EMCOOLS



EASY
EFFICIENT
PATIENT
COOLING



MECHANICAL CHEST COMPRESSION DEVICE



PREHOSPITAL COOLING

EMCOOLS FLEX.PAD

HYPOCARBON[®]


EMCOOLS



EASY
EFFICIENT
PATIENT
COOLING



- * Cooling rates of up to 3.3°C/h
- * Biocompatible material (skin- and environmentally friendly, non-toxic)




EMCOOLS FLEX.PAD

CONTINUOUS COOLING CHAIN

EMCOOLS




EASY
EFFICIENT
PATIENT
COOLING



EMCOOLS FLEX.PAD

PREPARATION


EMCOOLS



1

Preparation


- * Remove transportation- and storage box
- * Freeze Flex.Pad™ at -8 to -11°C horizontally
- * **Note:** This temperature range reflects a pre-cooling phase of 48h (during first-time freezing)



2

Operational readiness

- * Ready to use pads are indicated by a blue color indicator label
- * If the pad is too cold the label turns black



3

Application

- * Take out ready to use Flex.Pad™ (blue color indicator label)
- * Open packaging, remove protective foil and apply immediately on **dry skin**
- * **Note:** For activating the adhesive press down evenly for 3-5 seconds

PRE-HOSPITAL START WITH FLEX.PAD

INITIAL SITUATION - OUT-OF-HOSPITAL CARDIAC ARREST

EMCOOLS



EASY
EFFICIENT
PATIENT
COOLING



PRE-HOSPITAL TREATMENT (1/30/2015)

- * 16:25 69 year old male, loses consciousness in public bus
- * 16:26 Emergency call & start bystander basic life support
- * 16:33 Arrival of Vienna Ambulance Service & advanced cardiac life support
 - Ventricular fibrillation as initial rhythm
 - 6 mg epinephrine , 300 mg amiodarone
 - 6 defibrillations
- * 16:46 → ROSC
- * Application of Flex.Pad at an initial temperature of 35.3°C



PRE-HOSPITAL START WITH FLEX.PAD

EMCOOLS



1. Sedation, analgesia and paralysis

- * Midazolam 10 mg
- * Fentanyl 0.1 mg
- * Atracurium 25 mg
- * Bolus administered every 30 min

2. Continuous temperature monitoring

Esophageal temperature probe (inserted by a tubus)

3. Can be used on top of defibrillator pads



PRE-HOSPITAL START WITH FLEX.PAD

APPLICATION

EMCOOLS



EASY
EFFICIENT
PATIENT
COOLING




Flex.Pads are applied by the first responders

- * 1 pad on the chest (not on female breast tissue)
 - * 2 pad on the back
 - * 1 pad on the abdomen
 - * 1 pad per thigh
- The application procedure takes about 8 minutes

Rule-of-thumb

- * 1 Flex.Pad per 10 kg body weight



Notes

- * Don't apply on face, toes, fingers, genital region, female breast tissue or pregnant
- * Don't apply in case of skin diseases, inflammation, burns or any other skin injuries
- * During treatment skin temperature does not drop below 4 to 8°C

PRE-HOSPITAL COOLING

CLINICAL DATA

EMCOOLS




15 patients after cardiac arrest

Time to start of cooling: **12 min**

Time to target temperature: **70 min**

Δ Tesophagus ROSC-admission: **1.2° C**

Feasible and safe
No skin lesions



Out-of-hospital surface cooling to induce mild hypothermia in human cardiac arrest.
Uray et al. Resuscitation 2008.


PRE-HOSPITAL START WITH FLEX.PAD

INITIAL SITUATION - OUT-OF-HOSPITAL CARDIAC ARREST

EMCOOLS



EASY
EFFICIENT
PATIENT
COOLING



PRE-HOSPITAL START WITH FLEX.PAD

IN-HOSPITAL CONTINUATION / CATH LAB

EMCOOLS



CATH LAB (1/30/2015)

- * 17:16 Arrival at Cath Lab (Cardiology, Hanusch Hospital)
- * PTCA and two stents into circumflex artery

→ All Flex.Pads are to be removed at 34°C

→ Temperature on arrival: 33.8°C → Initial phase already completed

IN-HOSPITAL START WITH FLEX.PAD

CATH LAB – ANGIOGRAPHY DURING COOLING

EMCOOLS



→ Pads are radiolucent

* Good data on safety and feasibility of cooling during angiography, X-Ray, CT, and MRI

1



X-Ray

2

CT

3

MRI



PRE-HOSPITAL START WITH FLEX.PAD

IN-HOSPITAL CONTINUATION / CARDIAC CARE UNIT

EMCOOLS



CARDIAC CARE UNIT (1/30/2015)

* 18:05 Arrival at the CCU

→ Maintenance and rewarming

Therapeutic Hypothermia at $33^{\circ}\pm 1^{\circ}\text{C}$ for 24 hours (1/30/2015 – 1/31/2015)

* Continuous temperature monitoring by esophageal and bladder temperature probes (used simultaneously)

* Continuous administration of analgo-sedation and paralytics (standard medication)

PRE-HOSPITAL START WITH FLEX.PAD

IN-HOSPITAL CONTINUATION / CARDIAC CARE UNIT

EMCOOLS



EASY
EFFICIENT
PATIENT
COOLING



Maintenance at 33-34°C

- * 2 Flex.Pad reapplied twice to maintain the targeted temperature range

Passive Rewarming to 37°C

- * Patient only covered by a regular blanket
- * Target is rewarming to 37°C in about 8 hours, at a rewarming rate of 0.4°C/h
- * At 37.0°C all medication is stopped


PRE-HOSPITAL START WITH FLEX.PAD

IN-HOSPITAL CONTINUATION / CARDIAC CARE UNIT

EMCOOLS



Maintenance and Rewarming (1/31/2015)



PRE-HOSPITAL START WITH FLEX.PAD

IN-HOSPITAL CONTINUATION / CARDIAC CARE UNIT

EMCOOLS



EASY
EFFICIENT
PATIENT
COOLING



CCU (1/31/2015 – 2/02/2015)

Normothermia

* 37.0°C for day two and three (until 72 hours after cardiac arrest)

Fever control

* Paracetamol, Novalgin

* EMCOOLS Flex.Pad™ or other form of active temperature control


PRE-HOSPITAL START WITH FLEX.PAD

IN-HOSPITAL CONTINUATION / CARDIAC CARE UNIT

EMCOOLS



EASY
EFFICIENT
PATIENT
COOLING



PRE-HOSPITAL START WITH FLEX.PAD

EXAMINATIONS & FURTHER TREATMENT OF SEPSIS / PNEUMONIA

EMCOOLS




CCU 2/04/2015

Chest x-ray

- * Unremarkable, but putrid secretions from suctioning
- * Started Amoxi/Clav, changed to Moxifloxacin, finally Meropenem/Linezolid
- * Tracheal secretions
- * Pseudomonas aeruginosa,
- * Citrobacter freundii

→ **Pneumonia and sepsis**



THERAPEUTIC HYPOTHERMIA

SIDE EFFECTS - RISK OF INFECTION

EMCOOLS



Therapeutic Hypothermia and the Risk of Infection:

- * A Systematic Review and Meta-Analysis; Geurts et al, CCM 2013

Objectives

- * Systematic review and meta-analysis of randomized trials to examine the risk of infections in patients treated with hypothermia

Results

- * 23 studies, 2820 patients, 1396 cooled patients
- * Cooling for several hours-days (prodecures, cardiac arrest, TBI, stroke)

THERAPEUTIC HYPOTHERMIA

SIDE EFFECTS - RISK OF INFECTION

EMCOOLS



EASY
EFFICIENT
PATIENT
COOLING




Findings

All infections: no increased risk

Risk of pneumonia ↑
(risk ratio 1.44 [95% CI, 1.10–1.90])

Risk of sepsis ↑
(risk ratio 1.80 [95% CI, 1.04–3.10])

Overall prevalence of sepsis was low.



THERAPEUTIC HYPOTHERMIA

SIDE EFFECTS - RISK OF INFECTION

EMCOOLS



EASY
EFFICIENT
PATIENT
COOLING



Conclusio

- * Lack of definition of infections
- * Assessment of infections not blinded
- * Only one cardiac arrest study included
- * Take-home-message: high alertness towards signs of infection

THERAPEUTIC HYPOTHERMIA

PROPHYLACTIC ANTIBIOTICS?

EMCOOLS



Prophylactic antibiotics are associated with a lower incidence of pneumonia in cardiac arrest survivors treated with targeted temperature management[☆]

David J. Gagnon^{a,e,*}, Niklas Nielsen^{b,e}, Gilles L. Fraser^{a,c,e}, Richard R. Riker^{c,d,e}, John Dziodzio^{c,e}, Kjetil Sunde^{e,f}, Jan Hovdenes^{e,g}, Pascal Stammet^{e,h}, Hans Friberg^{e,i}, Sten Rubertsson^{e,j}, Michael Wanscher^{e,k}, David B. Seder^{c,d,e}

Univariate analysis of outcomes related to prophylactic antibiotics.

Variable	All Patients (n = 1240)	PRO (n = 416)	No-PRO (n = 824)	p ^a
Serious infections, n (%) ^b	(n = 1206)	(n = 414)	(n = 792)	
Any serious infection	553 (45.9%)	67 (16.2%)	486 (61.4%)	<0.001
Pneumonia	487 (40.4%)	52 (12.6%)	435 (54.9%)	<0.001
Sepsis	50 (4.1%)	5 (1.2%)	45 (5.7%)	<0.001
Other	57 (4.7%)	13 (3.1%)	44 (5.6%)	0.06
None	653 (54.1%)	347 (83.8%)	306 (38.6%)	<0.001

Conclusio

Prophylactic antibiotics were associated with a reduced incidence of pneumonia
Functional outcome was similar

PRE-HOSPITAL START WITH FLEX.PAD

NEUROLOGIC RECOVERY

EMCOOLS



CCU (02/03/2015)

- * Patient regained consciousness
- * Patient is extubated on 02/06/2015

Time on CCU: 01/30-02/11/2015

Regular Ward: 02/11-02/20/2015

→ **No neurologic sequelae (OPC 1, CPC 1)**

THERAPEUTIC HYPOTHERMIA

SEDATION

EMCOOLS



Sedation Standard - Department for Emergency Medicine, Medical Univ. Vienna

- * Midazolam 0,125 mg/kg/h
- * Fentanyl 0,002 mg/kg/h
- * Esmeron 0,25 mg/kg/h

Hypothermia and Drugs

- * Affects drug metabolism
- * Affects drug/receptor interaction
- * During hypothermia: high blood levels, reduced effect
- * During rewarming: toxicity develops

THERAPEUTIC HYPOTHERMIA

SEDATION

EMCOOLS



Comparison of midazolam/ fentanyl versus propofol/remifentanyl

- * Randomized clinical study comparing, 50 Patients
- * Lower time to offset for propofol/remifentanyl 13.2 (2.3–24) vs. 36.8 (28.5–45.1)
- * Norepinephrine infusion needed twice as often
- * Same outcome

* Free full text review:

* **Zhou and Poloyac**

**The effect of therapeutic hypothermia on drug metabolism and drug response:
cellular mechanisms to organ function**

* Expert Opin Drug Metab Toxicol. 2011 July ; 7(7): 803–816.

THERAPEUTIC HYPOTHERMIA

NEUROMUSCULAR BLOCKADE

EMCOOLS




Post hoc analysis of a prospective observational study:

- * Better survival
- * improved lactate clearance

Some problems

- * longer collapse to ROSC interval in the No-MB-group




THERAPEUTIC HYPOTHERMIA

EFFECT ON OUTCOME – VF PATIENTS, HACA-STUDY

EMCOOLS



First randomized trials (2002)




THERAPEUTIC HYPOTHERMIA

HACA-STUDY

EMCOOLS



EASY
EFFICIENT
PATIENT
COOLING



33°C vs. 36°C

How cool is cool enough?


EMCOOLS



EASY
EFFICIENT
PATIENT
COOLING




The NEW ENGLAND JOURNAL of MEDICINE




ORIGINAL ARTICLE

Targeted Temperature Management at 33°C versus 36°C after Cardiac Arrest



33°C vs. 36°C

How close do you want your brain to the heat?



33°C vs. 36°C

How close do you want your brain to be to the heat?

EMCOOLS




Figure 1. Mean temperature recordings for all patients. *Horizontal bars* mark the target temperature range of 32–34°C. Time 0 represents cooling initiation (n = 32).



ONGOING STUDIES IN:

- **Cooling during CPR**
- **Stroke**
- **MI**
- **Traumatic brain injury**
- **Spinal chord injury**
- **Hepatic encephalopathy**

Heat stroke



Thank you!